LASERPCINT

PcPlug-U

USB Communication Interface

(rev.03 - 11-07-2019)

Table of Contents:

1. In	stall	3
2. Se	et Up the Communication:	3
	ommands and Answers Formatting:	
3.1	Command format	
3.2	Answer format	
3.3	Error Message	
4. Pr	oduct Series and associated Command Tables	
	ommands and Answer Tables for Product Series #1	
5.1	Information commands	
5.2	Measurement setup commands	
5.3	Gain and Full Scale Commands	
5.4	Wavelength setup commands	
5.5	Measurement acquisition commands	
6. C	ommands and Answer Tables for Product Series #2	
6.1	Information commands	
6.2	Measurement setup commands	
6.3	Gain and Full Scale Commands	
6.4	Wavelength setup commands	12
6.5	Measurement acquisition commands	13
7. C	ommands and Answer Tables for Product Series #3	14
7.1	Information commands	14
7.2	Measurement setup commands	15
7.3	KEFUN=0-1-2-3-4 Gain and Full Scale Commands	
7.4	Wavelength setup commands	16
7.5	Measurement acquisition commands	16
8. Ex	camples and Notes	18
8.1	EXAMPLE 1 – quick communication test	18
8.2	EXAMPLE 2 – opening the communication using FTDXX functions	18
8.3	EXAMPLE 3 – command sequence to perform an energy measure	
8.4	EXAMPLE 4 – command sequence to perform a power measure	20
9. A :	nnex 1 - FTD2XX.DLL Dynamic Library	21
9.1	D2XX Driver Architecture	
9.2	Useful links for FTDI drivers and libraries:	22
9.3	DLL Functions	22

1. Install

Connect the **PcPlug-U** electronics through the USB port to the host PC.

Follow the prompt to install the USB drivers. They will be automatically installed through Windows update automatic download, or manually installed by using the driver installers available on the provided USB flash drive: inside the USB flash folder

X:\DRIVER\
you can find the driver installers for both 64 and 32 bit PCs:
\dpinst_amd64.exe
\dpinst_x86.exe

2. SET UP THE COMMUNICATION:

The **PcPlug-U** employs an FTDI FT232RQ chip that is a USB to serial UART interface.

This means that, for the host PC and for the user, this device will appear like a true COM port.

Example:

The COM number associated to this device can be known using WindowsTM "Device Manager" panel. And any serial terminal software (like WindowsTM HyperTerminal) can be used to open and test the communication with this device. See Chapter 8.1 for the example

These are the communication parameters that must be used:

Baud Rate:	For product Series #2 and #3: 38400 bps
	For product Series #1: 9600 kbps
	(See chapter 4 for part number lists)
Parity:	no
Data Bits:	8
Stop Bits:	1

When writing your Code, include the FTD2XX.DLL Dynamic Library for Windows in order to write your application. This library provides useful functions to recognize and set up the communication with the *PcPlug-U* device. Please see the FTD2XX.DLL complete manual on FTDI website, or refer to the extract we included in Annex 1 at the end of this paper.

3. COMMANDS AND ANSWERS FORMATTING:

When the *PcPlug-U* receives a valid input command, it confirms to the host device that the command has been received and return an answer. The Commands and Answers are ASCII formatted.

3.1 Command format

The format of a valid command is: *COMMANDNAME:

where:

"*"	Start of command	
"COMMANDNAME"	The command instruction is a sequence of ASCII characters. It must be in capitals.	
	Each command is described in the following pages.	
","	End of command	

3.2 Answer format

When the **PcPlug-U** receives a valid Input Command it replies sending a message through the USB interface.

The format of an answer is:

#ANSWER;

where:

"#"	Start of answer	
"ANSWER"	there are three kind of answer:	
	String: ASCII character sequence	
	• Int: integer number, numerical sequence (in ASCII code)	
	• Float: floating point number, numerical sequence plus decimal point (in ASCII	
	code)	
","	End of answer	

Please note that maximum response time from **PcPlug-U** is ~50msec. It's therefore recommended to set a delay of 50ms between write and read function.

3.3 Error Message

When the **PcPlug-U** receives an invalid command or a command affected by communication errors the answer is "??;"

where:

	"??"	USB communication error	
Ī	"," ,	End of answer	

Common reasons of error message are:

Input command not started with * character

Input command not in capitals

Input command does not correspond with the command list

4. PRODUCT SERIES AND ASSOCIATED COMMAND TABLES.

This manual includes different communication tables for different product series. Each product series has its specific commands and setting.

Please use the following table to find the correct commands and settings tables for your device.

Series #1 –		Seri	es #2 –	Series #3 –
OEM Series:		Thermopile Series:		BLINK Series:
		W-1500-D40-HPB-U W-1500-D40-SHC-U W-3000-D55-HPB-U W-3000-D55-SHC-U W-6000-D55-SHC-U		
For these devices refer to: Chapter 5.			evices refer to: pter 6.	For these devices refer to: Chapter 7.

5. COMMANDS AND ANSWER TABLES FOR PRODUCT SERIES #1

This product series includes these products:

Series #1 – OEM Series			
CSA-2-D12-BBF-U	Fit-50-H-U		
CSA-2-D12-HPB-U	Fit-200-H-U		
CSA-5-D12-BBF-U	Fit-500-H-U		
CSA-20-D20-BBF-U	Fit-3000-H-U		
CSA-20-D20-HPB-U	Fit-6000-H-U		
CSW-50-D20-BBF-U	Fit-IPL-R-H-U		
CSW-50-D20-HPB-U			
CSW-50-D25-BBF-U			
CSW-50-D25-HPB-U			
CSW-200-D20-HPB-U			
CSW-200-D30-HPB-U			

5.1 Information commands

These commands are used to get ID information about the **PcPlug-U** and the sensor. This info may be useful when asking Laserpoint for support.

Command	Answer (example)	Description
HEADN	"H" + String 8 char	Displays the Sensor Head model name (shortened)
SERNU	"S" + Int 6 digit	Displays the Sensor Head serial number
FHV	"H" + 2 char + "F" + 4	Displays the PcPlug-U Electronics Hardware and Firmware
	char	version
KEFUN	"K" + Int 2 digit	This 2 digit code number identifies the sensor type and the
		available functionalities.
		For this series of devices, the answer will be one among the
		highlighted:
		00 = OEM Thermopile sensor - Power
		01 = OEM Thermopile sensor - Fit Mode
		02 = OEM Thermopile sensor - Energy
		03 = OEM Thermopile sensor - Power + Energy
		04 = OEM Thermopile sensor - Fit Mode + Energy
		05 = Thermopile sensor — Power (see chapter 6)
		06 = Thermopile sensor - Power + Energy (see chapter 6)
		07 = Thermopile sensor — Fit mode (see chapter 6)
		08 = Thermopile sensor - Fit mode + Energy (see chapter 6)
		09 = Photodiode sensor
		10 = NA
		11 = NA
		12 = Blink Series Sensor - Power (see chapter 7)
		13 = Blink Series Sensor - Power + Energy (see chapter 7)

5.2 Measurement setup commands

These commands are used for the initial set up: operation mode selection (Power, Energy, others) and a reset of zero.

Command	Answer (example)	Description
POWER	"ok" or	Set PcPlug-U in Power Meter mode (if available)
	"NA" (if not available)	
ENERGY	"ok" or	Set PcPlug-U in Energy operation mode (if available)
	"NA" (if not available)	
ZERO	"ok"	Perform a Zero.
		This action will take about 3 seconds, and will reset the zero
		value of the sensor.
		Please make sure that this action is performed only when the
		sensor is not hit by laser or any other thermal source.
FAST	"FAST"	This is the default setting.
		It enables the acceleration algorithm, granting a faster
		response time.
SLOW	"SLOW"	Disables the acceleration algorithm.
		The response time may be dramatically lowered, but also noise
		(from laser source, or from cooling, or from environment) can
		be mitigated.
FASTSLOW	"FAST" or	Returns the Fast/Slow current setting
	"SLOW"	

5.3 Gain and Full Scale Commands

These commands allow selection of electronic amplifier gain (or in other words the selection of a Full Scale Range).

Command	Answer (example)	Description
SETX10	"ok"	Set the 0 th electronic amplifier gain (gain x1) (bigger full scale)
	"NA" (if not available)	
SETX1 1	"ok"	Set the 1 st electronic amplifier gain (gain x10) (smaller full scale)
	"NA" (if not available)	
X1D	Int 1 digit, 0 or 1	Displays the currently selected electronic amplifier gain:
		0: x1 gain
		1: x10 gain

5.4 Wavelength setup commands

Command	Answer (example)	Description
LAMBDA	"LAMBDA" + Int 1 digit	Displays the currently selected wavelength number.
	LAMBDA3	From number 1 to 5.
SETLAM + Int 1 digit	"ok"	Select the wavelength number.
SETLAM2		From number 1 to 5.
NOML1	String 3 char	Displays wavelength 1 Label
	CO2	
NOML2	String 3 char	Displays wavelength 2 Label
	YAG	
NOML3	String 3 char	Displays wavelength 3 Label
110107	LDS	
NOML4	String 3 char	Displays wavelength 4 Label
NOVE 5	VIS	Di la destala
NOML5	String 3 char EXC	Displays wavelength 5 Label
CFWL1	Float 2int.3dec	Displays the spectral correction coefficient of wavelength 1.
CFWLI	00.000	If the value is 00.000 it means that this wavelength is not
	00.000	available.
CFWL2	Float 2int.3dec	Displays the spectral correction coefficient of wavelength 2.
	01.000	If the value is 00.000 it means that this wavelength is not
		available.
CFWL3	Float 2int.3dec	Displays the spectral correction coefficient of wavelength 3.
	00.950	If the value is 00.000 it means that this wavelength is not
		available.
CFWL4	Float 2int.3dec	Displays the spectral correction coefficient of wavelength 4.
	00.990	If the value is 00.000 it means that this wavelength is not
		available.
CFWL5	Float 2int.3dec	Displays the spectral correction coefficient of wavelength 5.
	00.000	If the value is 00.000 it means that this wavelength is not
		available.

<u>NOTE:</u> in order to choose the correct wavelength is recommended to use all the NOMLx and CFWLx commands, so you can know which wavelengths are available (if CFWL answer \neq 0) or not available (if CFWL answer = 0).

5.5 <u>Measurement acquisition commands</u>

Command	Answer (example)	Description
OUTPM	Float	Displays measured power (W).
	4.325	The answer is formatted according to VISCA command.
VISCA	Int 1 digit	This command is used to know the measured number format
		and Unit of measure:
		0 = unit of measure W (or J) - no decimal number
		1 = unit of measure W (or J) – one decimal number
		2 = unit of measure W (or J) - two decimal number
		3 = unit of measure mW (or mJ) – no decimal number
		4 = unit of measure mW (or mJ) – one decimal number
		5 = unit of measure mW (or mJ) – two decimal number
		6 = unit of measure W (or J) – no decimal number (steps of 5W
		or 10W depending on sensor and gain)
STATUS	Int 3 digits	Displays status byte.
	114	Notice that this 3 digit integer must be converted into binary.
		bit 0: arm/zeroing done; (1) yes, (0) no
		bit 1: measure running; (1) yes, (0) no
		bit 2: Head connected; (1) yes, (0) no
		bit 3: cool alarm running; (1) yes, (0) no
		bit 4: wait before start a new measure; (1) yes
		bit 5: not used; default value (0)
		bit 6: overflow alarm; (1) yes, (0) no
		bit 7: thermistor connected; (1) yes, (0) no
TERMI	Int 1 digit	Thermistor availability: (1) yes, (0) no
		If thermistor is available the temperature
TEMP	Int 3 digit	Displays Head temperature x 10 (°C)
	255	

6. COMMANDS AND ANSWER TABLES FOR PRODUCT SERIES #2

This product series includes these products:

Series #2 - Thermopile Series			
A-02-D12-BBF-U	A-40/200-D40-HPB-U	W-200-D40-HPB-U	10-BB-D25-U
A-2-D12-BBF-U	A-40/200-D60-HPB-U	W-200-D40-SHC-U	10-UVA-D25-U
A-2-D12-HPB-U	A-200-D25-HPB-U	W-500-D70-SHC-U	10-UVC-D25-U
A-5-D12-BBF-U	A-200-D25-SHC-U	W-600-D30-HPB-U	20-BB-D40-U
A-10-D12-HPB-U	A-200-D40-HPB-U	W-600-D30-SHC-U	20-UVA-D40-U
A-10-D20-BBF-U	A-200-D40-SHC-U	W-1500-D40-HPB-U	20-UVC-D40-U
A-10-D20-HPB-U	A-200-D60-HPB-U	W-1500-D40-SHC-U	A-30-D12-SHC-L-U
A-30-D25-HPB-U	A-200-D60-SHC-U	W-3000-D55-HPB-U	10-BB-D12-L-U
A-40-D25-BBF-U	A-300-D60-HPB-U	W-3000-D55-SHC-U	A-10-D12-DIF-U
A-40-D25-HPB-U	A-600-D40-HPB-U	W-6000-D55-SHC-U	A-30-D18-DIF-U
A-40-D40-HPB-U	A-600-D60-SHC-U		A-40-D33-DIF-U
A-40/200-D25-HPB-U	A-1200-D60-SHC-U		

6.1 Information commands

These commands are used to get ID information about the **PcPlug-U** and the sensor. This info may be useful when asking Laserpoint for support.

Command	Answer (example)	Description
HEADN	"H" + String 8 char	Displays the Sensor Head model name (shortened)
SERNU	"S" + Int 6 digit	Displays the Sensor Head serial number
FHV	"H" + 2 char + "F" + 4	Displays the PcPlug-U Electronics Hardware and Firmware
	char	version
KEFUN	"K" + Int 2 digit	This 2 digit code number identifies the sensor type and the
		available functionalities.
		For this series of devices, the answer will be one among the
		highlighted:
		00 = OEM Thermopile sensor - Power (see chapter 5)
		01 = OEM Thermopile sensor - Fit Mode (see chapter 5)
		02 = OEM Thermopile sensor - Energy (see chapter 5)
		03 = OEM Thermopile sensor - Power + Energy (see chapter 5)
		04 = OEM Thermopile sensor - Fit Mode + Energy (see chapter 5)
		05 = Thermopile sensor – Power
		06 = Thermopile sensor - Power + Energy
		07 = Thermopile sensor – Fit mode
		08 = Thermopile sensor - Fit mode + Energy
		09 = Photodiode sensor
		10 = NA
		11 = NA
		12 = Blink Series Sensor - Power (see chapter 7)
		13 = Blink Series Sensor - Power + Energy (see chapter 7)

6.2 Measurement setup commands

These commands are used for the initial set up: operation mode selection (Power, Energy, others) and a reset of zero.

Command	Answer (example)	Description
POWER	"ok" or	Set PcPlug-U in Power Meter mode (if available)
	"NA" (if not available)	
ENERGY	"ok" or	Set PcPlug-U in Energy operation mode (if available)
	"NA" (if not available)	
ZERO	"ok"	Perform a Zero.
		This action will take about 3 seconds, and will reset the zero
		value of the sensor.
		Please make sure that this action is performed only when the
		sensor is not hit by laser or any other thermal source.
FAST	"FAST"	This is the default setting.
		It enables the acceleration algorithm, granting a faster
		response time.
SLOW	"SLOW"	Disables the acceleration algorithm. The response time may be
		dramatically lowered, but also noise (from laser source, or from
		cooling, or from environment) can be mitigated.
FASTSLOW	"FAST" or "SLOW"	Returns the Fast/Slow current setting

6.3 Gain and Full Scale Commands

These commands allow selection of electronic amplifier gain (or in other words the selection of a Full Scale Range).

Command	Answer (example)	Description
SETX1 0	"ok"	Set the 0 th electronic amplifier gain (biggest full scale)
	Or "NA" (if not available)	
SETX1 1	"ok"	Set the 1 st electronic amplifier gain
SETX1 2	"ok" "NA" (if not available)	Set the 2 nd electronic amplifier gain (smallest full scale)
SETX1 3	"ok" "NA" (if not available)	Set the "automatic" selection of the electronic amplifier gain
X1D	Int 1 digit, from 0 to 5	Displays the currently selected electronic gain set up:
		0: x1 gain
		1: x10 gain
		2: x100 gain
		3: automatic gain, current x1 gain
		4: automatic gain, current x10 gain
		5: automatic gain, current x100 gain
FSWX1 0	Number + _ + Unit of measure 10.0000_W	Command used to know the Power Full Scale of the 0 th gain. NOTE: this command and the following "FS*X1*" commands
	_	are also useful to get the formatting and the Unit of measure of
		the values received when using the acquisition commands (see "Measurement acquisition commands")
FSWX1 1	Number + _ + Unit of measure 5.0000_W	Command used to know the Power Full Scale of the 1st gain.
FSWX1 2	Number + _ + Unit of measure 1000.00_mW	Command used to know the Power Full Scale of the 2 nd gain.
FSJX1 0	Number + _ + Unit of measure NA	Command used to know the Energy Full Scale of the 0 th gain.
FSJX1 1	Number + _ + Unit of measure 10.0000_J	Command used to know the Energy Full Scale of the 1st gain.
FSJX1 2	Number + _ + Unit of measure 1000.00_mJ	Command used to know the Energy Full Scale of the 2 nd gain.

6.4 Wavelength setup commands

Command	Answer (example)	Description
LAMBDA	LAMBDA + Int 5 digit	Displays the currently selected wavelength (nm)
	LAMBDA01064	
SETLAM + Int 5 digit	"LAMBDA" + Int 5 digit	Command used to set the desired wavelength (nm).
SETLAM00970	LAMBDA00970	
RANGEWL	"RWL_" + Int5digit +	Displays the minimum and maximum value of the wavelength
	"_to_" + Int5digit	range (nm). Any wavelength among this range can be
	RWL_00200_to_01100	selected.
		(see SETLAM command)
SINGLEWL	"SWL_" + Int5digit + "_"	Displays a variable list of wavelengths (nm).
	+ Int5digit + "_" +	These are discrete wavelength and only the specific value
	SWL_1550_2940_10600	values are wavelength that can be selected (see SETLAM
		command).

6.5 Measurement acquisition commands

Command	Answer (example)	Description
OUTPM	Float	Displays measured power value (or energy).
	4.325	The number format varies depending on many parameters. To
		know the number format of each scale, please use the "FS*X1
		*" commands.
		This is a "one command – one answer" command: each time
		this command is sent, one one value is answered.
		This command is used to request a few samples per seconds
		(max 5-8 requests). For higher sampling it's recommended to
		use the OUTPTS command.
STATUS	"Y" + Int 5 digits	Displays the status byte.
	Y00003	Notice that this 5 digit integer must be converted into binary.
		Bit 0: Head connected: (1) yes, (0) no
		Bit 1: thermistor connected: (1) yes, (0) no
		Bit 2: not used
		Bit 3: cool warning (1)
		Bit 4: battery: connected to AC (1)
		Bit 5: battery: charge in progress (1)
		Bit 6: overload warning (1)
		Bit 7: overflow warning (1)
		Bit 8: status "ready", for Fit/Energy mode (1)
		Bit 9: status "triggered", for Fit/Energy mode (1)
		Bit 10: status "wait", for Fit modes (1)
		Bit 11: not used
		Bit 12: overflow ADC gain G=x1 (1)
		Bit 13: overflow ADC gain G=x10 (1)
		Bit 14: overflow ADC gain G=x100 (1)
		Bit 15: not used
TERM	"T" + Int 1 digit	Thermistor availability: (1) yes, (0) no
TEMP	"t" + Int 3 digit	Displays Head temperature x 10 (°C)
	t258	
OUTPTS	Integer + "_" + Integer +	This command is used to get a continuous stream of measured
	"_" + Integer	values (formatted like OUTPM command) + status byte
	0.0994_00003_258	(formatted like STATUS) + sensor temperature (formatted like TEMP).
		It is used when power mode is selected and a continuous flow
		of measurement is required.
		The output is delivered at 8hz, meaning 125ms interval between
		two measures.

7. COMMANDS AND ANSWER TABLES FOR PRODUCT SERIES #3

This series includes these products:

Series #3 – BLINK Series			
BL-W-50W-16-K-U			
BL-A-30W-16-K-U			

Please notice that the Baud Rated for this product series is 38400 bps (see Chapter 2 for communication setup)

7.1 Information commands

These commands are used to get ID information about the **PcPlug-U** and the sensor. This info may be useful when asking Laserpoint for support.

Command	Answer (example)	Description
HEADN	"H" + String 8 char	Displays the Sensor Head model name (shortened)
SERNU	"S" + Int 6 digit	Displays the Sensor Head serial number
FHV	"H" + 2 char + "F" + 4	Displays the PcPlug-U Electronics Hardware and Firmware
	char	version
KEFUN	"K" + Int 2 digit	This 2 digit code number identifies the sensor type and the available functionalities. For this series of devices, the answer will be one among the
		highlighted:
		00 = OEM Thermopile sensor - Power (see chapter 5)
		01 = OEM Thermopile sensor - Fit Mode (see chapter 5)
		02 = OEM Thermopile sensor - Energy (see chapter 5)
		03 = OEM Thermopile sensor - Power + Energy (see chapter 5)
		04 = OEM Thermopile sensor - Fit Mode + Energy (see chapter 5)
		05 = Thermopile sensor — Power (see chapter 6)
		$\theta = \text{Thermopile sensor} - \text{Power} + \text{Energy}$ (see chapter 6)
		07 = Thermopile sensor — Fit mode (see chapter 6)
		08 = Thermopile sensor - Fit mode + Energy (see chapter 6)
		09 = Photodiode sensor
		10 = NA
		11 = NA
		12 = Blink Sensor - Power
		13 = Blink Sensor - Power + Energy

7.2 Measurement setup commands

These commands are used for the initial set up: operation mode selection (Power, Energy, others) and a reset of zero.

Command	Answer (example)	Description
POWER	"ok" or	Set PcPlug-U in Power Meter mode (if available)
	"NA" (if not available)	
ENERGY	"ok" or	Set PcPlug-U in Energy operation mode (if available)
	"NA" (if not available)	
ZERO	"Zok"	Perform a Zero.
		This action will take about 3 seconds, and will reset the zero
		value of the sensor.
		Please make sure that this action is performed only when the
		sensor is not hit by laser or any other thermal source.
FAST	"FAST"	This is the default setting.
		It enables the acceleration algorithm, granting a faster
		response time.
SLOW	"SLOW"	Disables the acceleration algorithm. The response time may be
		dramatically lowered, but also noise (from laser source, or from
		cooling, or from environment) can be mitigated.
FASTSLOW	"FAST" or "SLOW"	Returns the Fast/Slow current setting

7.3 KEFUN=0-1-2-3-4 Gain and Full Scale Commands

These commands allow selection of electronic amplifier gain (or in other words the selection of a Full Scale Range).

Command	Answer (example)	Description
SETX1 0	"ok"	Set the 0 th electronic amplifier gain (biggest full scale)
	"NA" (if not available)	
SETX1 1	"ok"	Set the 1st electronic amplifier gain
	or "NA" (if not available)	
SETX1 2	"ok"	Set the 2 nd electronic amplifier gain (smallest full scale)
	"NA" (if not available)	
SETX1 3	"ok"	Set the "automatic" selection of the electronic amplifier gain
	"NA" (if not available)	
X1D	Int 1 digit, from 0 to 5	Displays the currently selected electronic gain set up:
		0: x1 gain
		1: x10 gain
		2: x100 gain
		3: automatic gain, current x1 gain
		4: automatic gain, current x10 gain
		5: automatic gain, current x100 gain
FSWX1 0	Number + _ + Unit of measure	Command used to know the Power Full Scale of the 0 th gain.
	10.0000_W	NOTE: this command and the following "FS*X1*" commands
		are also useful to get the formatting and the Unit of measure of
		the values received when using the acquisition commands (see
		"Measurement acquisition commands")
FSWX1 1	Number + _ + Unit of measure	Command used to know the Power Full Scale of the 1st gain.
	5.0000_W	
FSWX1 2	Number + _ + Unit of measure 1000.00_mW	Command used to know the Power Full Scale of the 2 nd gain.
FSJX1 0	Number + _ + Unit of measure	Command used to know the Energy Full Scale of the 0 th gain.

	NA	
FSJX1 1	Number + _ + Unit of measure 10.0000_J	Command used to know the Energy Full Scale of the 1st gain.
FSJX1 2	Number + _ + Unit of measure 1000.00_mJ	Command used to know the Energy Full Scale of the 2 nd gain.

7.4 Wavelength setup commands

Command	Answer (example)	Description
LAMBDA	LAMBDA + Int 5 digit	Displays the currently selected wavelength (nm)
	LAMBDA01064	
SETLAM + Int 5 digit	"LAMBDA" + Int 5 digit	Command used to set the desired wavelength (nm).
SETLAM00970	LAMBDA00970	
RANGEWL	"RWL_" + Int5digit +	Displays the minimum and maximum value of the available
	"_to_" + Int5digit	wavelength range (nm). Any wavelength among this range
	RWL_00200_to_01100	can be selected.
		(see SETLAM command)
SINGLEWL	"SWL_" + Int5digit + "_"	Displays a variable list of available wavelengths (nm).
	+ Int5digit + "_" +	These are discrete wavelength and only the specific value
	SWL_1550_2940_10600	values are wavelength that can be selected
		(see SETLAM command).

7.5 <u>Measurement acquisition commands</u>

Command	Answer (example)	Description
OUTPM	Float	Displays measured power value (or energy).
	4.325	The number format varies depending on many parameters. To
		know the number format of each scale, please use the "FS*X1
		*" commands.
		This is a "one command – one answer" command: each time
		this command is sent, one one value is answered.
		This command is used to request a few samples per seconds
		(max 5-8 requests). For higher sampling it's recommended to
		use the OUTPTS command.
OUTPTS	16 x [Integer + "_"] + "s"	This command is used to activate a continuous stream of
	+ 5 digit Integer + "t" + 3	measured values. It is used when power mode is selected and a
	digit Integer + "c" + 2	continuous flow of measurement with high sampling rate is
	digit Integer	required.
	3.056_3.054_3.052_3.049_	
	3.047_3.045_3.043_3.041_	The answer is delivered 12 times per second in form a of a
	3.038_3.036_3.034_3.032_ 3.030_3.028_3.026_3.025	string.
	s00003t251c49	Each string contains:
	3000031231049	16 measured values (formatted like OUTPM command) +
		Status bytes (formatted like STATUS command) +
		Sensor temperature (formatted like TEMP) +
		2 digit counter (increasing from 00 to 99)
		NOTE: The total number of values outputted per second is $= 12$
		strings x 16 values = 192 samples
		NOTE2: the counter that can be used to check if strings are
		complete and/or if there is some string that is missing due to
		communication errors.

		To stop the data stream use the command *COMMAND:
COMMAND	"COMMAND"	This command is used to stop any data stream mode that is active.
STATUS	"Y" + Int 5 digits Y00003	Displays the status byte. Notice that this 5 digit integer must be converted into binary. Bit 0: Head connected: (1) yes, (0) no Bit 1: thermistor connected: (1) yes, (0) no Bit 2: not used Bit 3: cool warning (1) Bit 4: battery: connected to AC (1) Bit 5: battery: charge in progress (1) Bit 6: overload warning (1) Bit 7: overflow warning (1) Bit 8: status "ready", for Fit/Energy mode (1)
		Bit 9: status "triggered", for Fit/Energy mode (1) Bit 10: status "wait", for Fit modes (1) Bit 11: not used Bit 12: overflow ADC gain G=x1 (1) Bit 13: overflow ADC gain G=x10 (1) Bit 14: overflow ADC gain G=x100 (1) Bit 15: not used
TERM TEMP	"T" + Int 1 digit "t" + Int 3 digit t258	Thermistor availability: (1) yes, (0) no Displays Head temperature x 10 (°C)

8. EXAMPLES AND NOTES

8.1 EXAMPLE 1 – quick communication test

The simplest way to test the communication, between the PC and a **PcPlug-U**, is using a simple serial terminal software (like Hyperterminal, Putty, CoolTerm, and many others).

Once the **PcPluq-U** has been connected to the PC and the drivers have been installed:

- set the correct COM port in the serial terminal software (use the device manager to get the COM number)
- use these settings:

Baud Rate: 38400 bps
Parity: no
Data Bits: 8
Stop Bits: 1

Start the connection and send a simple command:

*SERNU: to get a 6 digit serial number as an answer or *OUTPM: to get a measured value as an answer.

8.2 EXAMPLE 2 – opening the communication using FTDXX functions

Here below is reported an example of the main steps necessary to start the communication with a **PcPlug-U** device by using the FTDXXX functions. The programming language for this example is VB.NET.

```
'Get the number of the connected FTDI devices:

FT_Status = FT_GetNumberOfDevices(NumDevicesConnected, vbNullChar, FT_LIST_NUMBER_ONLY)

Browse all the connected FTDI devices to find Laserpoint PcPlug-U device
```

'1 Get the device description

FT_Status = FT_GetDeviceString(i, Description, FT_LIST_BY_INDEX Or FT_OPEN_BY_DESCRIPTION)

'2 Shrink the description returned as 64 chars string to the correct number of chars.

'Ex: " PcPlug_V3 " --> " PcPlug_V3"

Description = Microsoft.VisualBasic.Left(Description, InStr(1, Description, vbNullChar) - 1)

NOTE: Description of Laserpoint *PcPlug-U* devices is "*PcPlug_V3*"

Get the serial number of Laserpoint device using the description

```
'Get serial number of device using index
FT_Status = FT_GetDeviceString(i, Serial, FT_LIST_BY_INDEX Or FT_OPEN_BY_SERIAL_NUMBER)

'Shrink the description from 64 chars to the correct number of chars
'Ex: "123456"

Serial = Microsoft.VisualBasic.Left(Serial, InStr(1, Serial, vbNullChar) - 1)
```

FT_Status = FT_OpenBySerialNumber(Serial, FT_OPEN_BY_SERIAL_NUMBER, COM_Handle)

```
'Setting communication parameters

FT_Status = FT_SetBaudRate(COM_Handle, 38400)

FT_Status = FT_SetDataCharacteristics(COM_Handle, FT_DATA_BITS_8, FT_STOP_BITS_1, FT_PARITY_NONE)

FT_Status = FT_SetFlowControl(COM_Handle, FT_FLOW_NONE, 0, 0)
```

Open communication with device identified by its serial number. Function will return a communication Handle which will be used for all the following communications

```
FT_Status = FT_SetTimeouts(COM_Handle, 500, 100)
FT_Status = FT_Purge(COM_Handle, FT_PURGE_RX Or FT_PURGE_TX)
```

FT_Status = FT_Write_String(COM_Handle, "*COMMAND:", 9, byteswritten)

' Set a delay of 50ms between write and read function to wait for device to reply

System.Threading.Thread.Sleep(50)

' Read the response from the device

FT_Status = FT_GetQueueStatus(COM_Handle, RXBytes)

ReadString = Space(RXBytes)

FT_Status = FT_Read_String(COM_Handle, ReadString, RXBytes, bytesread)

8.3 EXAMPLE 3 – command sequence to perform an energy measure

Here below is reported an example of **PcPlug-U** command sequence to perform an ENERGY measure using a sensor from series #1 (sensors with **KEFUN code 0-1-2-3-4**)

Please note that the answers to the command "*STATUS:" are just an example and may be different depending on the sensor employed.

Command	Answer	Comment	
ENERGY	Ok	Activating energy measurement mode	
SETX1 0	Ok	Activating the "standard" full scale range	
NOML2	"YAG"	Checking the "name" of the wavelength saved in the 2 nd position of memory.	
CFWL2	0.982	Checking the coefficient of this wavelength (It's important to check that is not = 0. The "0" value means that this wavelength is not activated)	
You n availd		ommands, for each number from 1 to 5, in order to know which wavelength are	
SETLAM2	ok	Selecting the 2 nd wavelength	
STATUS	132 (=1000010 <u>0</u>)	As expected, the first time the device is started the answer is a "no" (=0) on bit number 0. That's because the instrument has not been zeroed/armed yet	
ZERO	ok	Zeroing the PcPlug-U	
		the sensor is in steady state while performing the zero, and no laser or thermal sor during this 3-4 seconds operation.	
STATUS	133 (=1000010 <u>1</u>)	Now the bit number 0 is =1 (=yes, the instrument has been zeroed). The PcPlug-U is armed and ready to measure.	
Turn	the laser ON and sho	ot.	
NOTE	: Laser pulse maximu	m duration is 100-300 ms depending on sensor type.	
STATUS	134 (=100001 <u>1</u> 0)	Bit number 1 is =1. Meaning that the measure is running. (this status will lasts 4-5 seconds usually)	
STATUS	148 (=100 <u>1</u> 0100)	Bit number 5 is = 1. Meaning "Wait before start new measure". (this status may last from 4 to 40 seconds, depending on sensor type). You may use the command *HOFTE: to know approximately the duration of this status (seconds).	
No LA	NTION: SER pulses must hit t ew measure").	he sensor during this span of time (while bit 5 = 1 = meaning "wait before start	
OUTPM	1.65	Getting the energy measured (Joules) NOTE: this value is available from the moment the "measure running" stops, and it will be available until a new measure starts or a new zero is performed	

^{&#}x27;Read and write to serial port

STATUS	132 (=1000010 <u>0</u>)	The waiting time is finished, The PcPlug-U is not yet armed.
STATUS	133 (=1000010 <u>1</u>)	The PcPlug is ready to run a new measure

8.4 EXAMPLE 4 – command sequence to perform a power measure

This is an example of **PcPlug-U** command sequence to perform a POWER measure using a sensor **of KEFUN family** 5-6-7-8-9 or 12-13.

Please note that the answers to the command "*STATUS:" are just an example and may be different depending on the model of employed sensor.

Command	Answer	Comment
RANGEWL	RWL_00200_to_01100	Get info about available wavelengths.
		This answer means that any wavelength in the range 200nm – 1100 nm
		can be selected.
SINGLEWL	SWL_1550_2940	This answer means that 1550nm and 2940nm are available as discrete values (e.g. 1600nm is not selectable).
SETLAM01070	LAMBDA01070	The wavelength of 1070nm has been selected.
FSWX1 0	20.0000_W	Command used to know the maximum Power and the unit of measure of the 0^{th} electronic gain (this Full Scale).
FSWX1 1	5.0000_W	Command used to know the maximum Power and the unit of measure of the 1st electronic gain (this Full Scale).
SETX1 1	ok	Select the 5W Full Scale
STATUS	Y00003 =	This status shows that. Bit 0: Head connected: (1) yes
	000000000000011	Bit 1: thermistor connected: (1) yes no alarms or warnings are activated.
OUTPM	0.0027	The measured value is 0.0027 Watts. It has to be evaluated an instrument Zeroing. Note that the number of digits outnumbers the precision and noise of the instrument: it's therefore recommended to round up the number to match the display needs.
ZERO	Zok	Perform a zero reset
OUTPM	0.0006	The measured value is now 0.0006 Watts. Note that the number of digits outnumbers the precision and the noise of the instrument. Therefore: - after zeroing, the answer will not be 0.0000 - it's recommended to round up this value according to the displaying needs
A	After the laser has been turn	
OUTPM	2.4986	The measured value is now 2.499 Watts. (rounded)


9. ANNEX 1 – FTD2XX.DLL DYNAMIC LIBRARY

The FTD2XX.DLL Dynamic Library for Windows allows you to write your application.

The architecture of the FTD2XX.DLL drivers consists of a Windows WDM driver that communicates with the device via the Windows USB Stack and a DLL which interfaces the Application Software (written in VC++, C++ Builder, Delphi, VB etc.) to the WDM driver.

The FTD2XX.DLL interface provides a simple, easy to use, set of functions to access **PcPlug-U** control card.

9.1 <u>D2XX Driver Architecture</u>



9.2 Useful links for FTDI drivers and libraries:

For a complete list and description of FTDXXX functions, please download the "FTD2XX Programmer's Guide" at this link:

http://www.ftdichip.com/Support/Documents/ProgramGuides.htm

To download the right libraries for your Operative System / architecture please check this section of FTDI website: http://www.ftdichip.com/Drivers/D2XX.htm

For software examples with different programming languages please check this section: http://www.ftdichip.com/Support/SoftwareExamples/CodeExamples.htm

9.3 DLL Functions

Here below we report only an extract of the "FTD2XX Programmer's Guide" (for the full document see the link above), that collects the most important functions that may be used to identify and communicate with the **PcPlug-U**:

FT ListDevices

Description Gets information concerning the devices currently connected. This function can return such information as the number of devices connected, and device strings such as serial number and product description.

Syntax FT_STATUS FT_ListDevices (PVOID pvArg1, PVOID pvArg2, DWORD dwFlags)

Parameters

pvArg1 meaning depend on the dwFlags value (see note below)pvArg2 meaning depend on the dwFlags value (see note below)dwFlags Determines format of returned information (see note below)

Return Value FT_OK if successful, otherwise the return value is an FT error code

Note Remarks This function can be used in a number of ways to return different types of information.

In its simplest form, it can be used to return the number of devices currently connected. If FT_LIST_NUMBER_ONLY bit is set in dwFlags, the parameter pvArg1 is interpreted as a pointer to a DWORD location to store the number of devices currently connected.

It can be used to return device string information. If

FT_OPEN_BY_SERIAL_NUMBER bit is set in dwFlags, the serial number string will be returned from this function. If FT_OPEN_BY_DESCRIPTION bit is set in dwFlags, the product description string will be returned from this function. If neither of these bits is set, the serial number string will be returned by default. It can be used to return device string information for a single device. If FT_LIST_BY_INDEX bit is set in dwFlags, the parameter pvArg1 is interpreted as the index of the device, and the parameter pvArg2 is interpreted as a pointer to a buffer to contain the appropriate string. Indexes are zerobased, and the error code FT_DEVICE_NOT_FOUND is returned for an invalid index.

It can be used to return device string information for all connected devices. If FT_LIST_ALL bit is set in dwFlags, the parameter pvArg1 is interpreted as a pointer to an array of pointers to buffers to contain the appropriate strings, and the parameter pvArg2 is interpreted as a pointer to a DWORD location to store the number of devices currently connected. Note that, for pvArg1, the last entry in the array of pointers to buffers should be a NULL pointer so the array will contain one more location than the number of devices connected.

FT_Open

Description Opens the device and return a handle which will be used for subsequent accesses.

Syntax FT_STATUS FT_Open (int iDevice, FT_HANDLE *ftHandle)

Parameters

iDevice indicates the number of the device to be opened. Must be 0 if only one device is attached. For multiple devices 1, 2 etc. *ftHandle* Pointer to a variable of type FT_HANDLE where the handle will be stored. This handle must be used to access the device.

Return Value FT_OK if successful, otherwise the return value is an FT error code

Note Although this function can be used to open multiple devices by setting iDevice to 0, 1, 2 etc. there is no ability to open a specific device. To open named devices, use the function **FT_OpenEx**. With the **FT_OpenEx** function (not described in this user manual) it is possible to open a device also trough its *serial number* or trough its description. For further information, please contact **LASERPOINT.srl**.

FT_Close

Description Closes the communication with a open device.

Syntax FT_STATUS FT_Close (FT_HANDLE ftHandle)

Parametres

ftHandle pointer to the communication handle of the device to close.

Return Value FT_OK if successful, otherwise the return value is an FT error code

FT_Read

Description Reads a string from the device.

Syntax FT_STATUS **FT_Read** (FT_HANDLE *ftHandle*, LPVOID *lpBuffer*, DWORD *dwBytesToRead*, LPDWORD *lpdwBytesReturned*)

Parameters

ftHandle pointer to the communication handle of the device to read.
 IpBuffer pointer to the buffer that receives the data from the device.
 DwBytesToRead Number of bytes to be read from the device.
 IpdwBytesReturned Pointer to a variable of type DWORD which receives the number of bytes read from the device.

Return Value FT_OK if successful, FT_IO_ERROR otherwise.

Note FT_Read always returns the number of bytes read in IpdwBytesReturned. This function does not return until dwBytesToRead have been read into the buffer. The number of bytes in the receive queue can be determined by calling FT_GetStatus or FT_GetQueueStatus, and passed to FT_Read as dwBytesToRead so that the function reads the device and returns immediately. When a read timeout value has been specified in a previous call to FT_SetTimeouts, FT_Read returns when the timer expires or dwBytesToRead have been read, whichever occurs first. If the timeout occurred, FT Read reads available data into the buffer and returns FT OK. An application should use the function return value and IpdwBytesReturned when processing the buffer. If the return value is FT_OK, and lpdwBytesReturned is equal to dwBytesToRead then FT_Read has completed normally. If the return value is FT_OK, and IpdwBytesReturned is less then dwBytesToRead then a timeout has occurred, and the read has been partially completed. Note that if a timeout occurred and no data was read, the return value is still FT OK. A return value of FT IO ERROR suggests an error in the parameters of the function, or a fatal error like USB disconnect has occurred.

FT_Write

Description Writes a string to the device.

Syntax FT_STATUS **FT_Write** (FT_HANDLE *ftHandle*, LPVOID *lpBuffer*, DWORD *dwBytesToWrite*, LPDWORD *lpdwBytesWritten*)

Parameters

ftHandle pointer to the communication *handle* of the device to write. **IpBuffer** pointer to the buffer which contains the bytes to be written in the device

DwBytesToWrite number of bytes to write to the device.

IpdwBytesWritten pointer to a variable of type DWORD which receives the number of bytes written to the device

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FT_ResetDevice

Description Sends a Reset command to the device.

Syntax FT_STATUS FT_ResetDevice (FT_HANDLE ftHandle)

Parameters

ftHandle pointer to the communication handle of the device to reset .

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FT SetBaudRate

Description Sets the baudrate for the device.

Syntax FT_STATUS **FT_SetBaudRate** (FT_HANDLE *ftHandle*, DWORD *dwBaudRate*)

Parameters

FtHandle pointer to the communication *handle* of the device to set out. **dwBaudRate** value of the *baudrate* to set out.

Return Value FT OK if successful, otherwise the return value is an FT error code.

Note: PcPlug-U Baud Rate value is 38400.

FT_SetDataCharacteristics

Description Sets the data characteristics for the device.

Syntax FT_STATUS **FT_SetDataCharacteristics** (FT_HANDLE *ftHandle*, UCHAR *uWordLength*, UCHAR *uStopBits*, UCHAR *uParity*)

Parameters

ftHandle pointer to the communication <code>handle</code> of the device to set out . **uWordLength** number of <code>bits</code> per word. It must set as <code>FT_BITS_8</code> (in the case of 8 bit schosen) or as <code>FT_BITS_7</code> (in the case of 7 bits chosen). **uStopBits** number of stop <code>bits</code>. It must set as <code>FT_STOP_BITS_1</code> (when one stop bit is requested) or as <code>FT_STOP_BITS_2</code> (when two stop bits are requested). **uParity** number of parity <code>bits</code>. It must set as <code>FT_PARITY_NONE</code> (no parity bit) or as <code>FT_PARITY_ODD</code> (parity bit is odd) or as <code>FT_PARITY_EVEN</code> (parity bit is even) or as <code>FT_PARITY_MARK</code> (always high parity bit) or as <code>FT_PARITY_SPACE</code> (always low parity bit).

Return Value FT_OK if successful, otherwise the return value is an FT error code.

Note: for **PcPlug-U** the DataCharacteristics must be set as FT_DATA_BITS_8, FT_STOP_BITS_1, FT_PARITY_NONE

FT_SetFlowControl

Description Sets the flow control the chip serial communication of chip USB/RS232.

Syntax FT_STATUS **FT_SetDataCharacteristics** (FT_HANDLE *ftHandle*, USHORT *usFlowControl*, UCHAR *uXon*, UCHAR *uXoff*)

Parameters

FtHandle pointer to the communication *handle* of the device to set out. **usFlowControl** set the kind of flow control. It must be set as FT_FLOW_NONE (no flow control) or as FT_FLOW_RTS_CTS (hardware RTS/CTS flow control) or as FT_FLOW_DTR_DSR (hardware DTR/DSR flow control) or as FT_FLOW_XON_XOFF (software XON/XOFF flow control) **uXon** shows the character uses as Xon signal. It must be set only when the flow control is software XON/XOFF kind (otherwise, it must be set as zero).

uXoff shows the character uses as Xoff signal. It must be set only when the flow control is software XON/XOFF kind (otherwise, it must be set as zero).

Return Value FT_OK if successful, otherwise the return value is an FT error code.

Note: for PcPlug-U the FlowControl must be set as FT FLOW NONE

FT_SetDTR

Description Sets the Data Terminal Ready (DTR) control signal. (Data Terminal Ready).

Syntax FT_STATUS FT_SetDTR (FT_HANDLE ftHandle)

Parameters

ftHandle pointer to the communication handle of the DTR device to set out.

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FT_CIrDTR

Description This function clears the Data Terminal Ready (DTR) control signal (*Data Terminal Ready*).

Syntax FT_STATUS **FT_CIrDTR** (FT_HANDLE *ftHandle*)

Parameters

ftHandle pointer to the communication handle of the DTR device to set out.

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FT_SetRTS

Description Sets the Request To Send (RTS) control signal. (Request To Send).

Syntax FT_STATUS FT_SetDTR (FT_HANDLE ftHandle)

Parameters

ftHandle pointer to the communication handle of the RTS device to set out.

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FT_CIrRTS

Description Clears the Request To Send (RTS) control signal (Request To Send).

Syntax FT_STATUS FT_SetDTR (FT_HANDLE ftHandle)

Parameters

FtHandle pointer to the communication handle of the RTS device to set out.

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FT_SetTimeouts

Description Sets the read and write timeouts for the device.

Syntax FT_STATUS **FT_SetBaudRate** (FT_HANDLE *ftHandle*, DWORD *dwReadTimeout*, DWORD *dwWriteTimeout*)

Parameters

FtHandle pointer to the communication *handle* of the device to set out . **dwReadTimeout** value of the Read timeout, in milliseconds, to set out. **dwWriteTimeout** value of the Write timeout, in milliseconds, to set out.

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FT_GetQueueStatus

Description Shows the number of characters in the receive queue.

Syntax FT_STATUS **FT_GetQueueStatus** (FT_HANDLE *ftHandle*, LPDWORD *lpdwAmountlnRxQueue*)

Parameters

FtHandle pointer to the communication handle of the device to set out .

IpdwAmountInRxQueue Pointer to a variable of type DWORD which receives the number of characters in the receive queue.

Return Value FT OK if successful, otherwise the return value is an FT error code.

FT_GetStatus

Description Shows the device status including number of characters in the receive queue, number of characters in the transmit queue, and the current event status.

Syntax FT_STATUS **FT_GetStatus** (FT_HANDLE *ftHandle*, LPDWORD *lpdwAmountlnRxQueue*, LPDWORD *lpdwAmountlnTxQueue*, LPDWORD *lpdwEventstatus*)

Parameters

ftHandle pointer to the communication handle of the device to set out .

IpdwAmountInRxQueu Pointer to a variable of type DWORD which receives the number of characters in the receive queue.

LpdwAmountInTxQueue Pointer to a variable of type DWORD which receives the number of characters in the transmit queue.

IpdwEventstatus Pointer to a variable of type DWORD which receives the current state of the event status.

Return Value FT_OK if successful, otherwise the return value is an FT error code.

FTD2XX - Error codes

```
FT OK = 0
FT_INVALID_HANDLE = 1
FT_DEVICE_NOT_FOUND = 2
FT_DEVICE_NOT_OPENED = 3
FT_IO_ERROR = 4
FT_INSUFFICIENT_RESOURCES = 5
FT_INVALID_PARAMETER = 6
FT_INVALID_BAUD_RATE = 7
FT_DEVICE_NOT_OPENED_FOR_ERASE = 8
FT_DEVICE_NOT_OPENED_FOR_WRITE = 9
FT FAILED TO WRITE DEVICE = 10
FT_EEPROM_READ_FAILED = 11
FT_EEPROM_WRITE_FAILED = 12
FT EEPROM ERASE FAILED = 13
FT_EEPROM_NOT_PRESENT = 14
FT_EEPROM_NOT_PROGRAMMED = 15
FT_INVALID_ARGS = 16
```